

PACKIT

***T. foetus* Reagent Set**

For *Tritrichomonas foetus* Detection

User Manual

For Research Use Only

Manufacturer:

GeneReach Biotechnology Corporation

TEL: 886-4-24639869 FAX: 886-4-24638255

No. 19, Keyuan 2nd Rd., Central Taiwan Science Park, Taichung City 407, Taiwan

Web Site: www.genereach.com

2015/08 d

Content

INTENDED USE.....1

SUMMARY AND EXPLANATION.....1

PRINCIPLES OF THE PROCEDURE.....2

PRODUCT DESCRIPTION.....2

A. Materials Provided.....2

B. Materials and Equipment Required, but Not Provided.....3

C. Storage and Stability.....3

PRECAUTIONS.....3

LIMITATIONS.....4

SAMPLE TYPE.....5

OPERATION PROCEDURE.....5

A. Using **POCKIT™ *T. foetus*** Reagent Set.....5

B. Quick Guide.....7

DATA INTERPRETATION.....8

ANALYTICAL SENSITIVITY.....8

TROUBLESHOOTING.....9

REFERENCE.....11

INTENDED USE

POCKIT™ *T. foetus* Reagent Set uses insulated isothermal polymerase chain reaction (iiPCR) technology (Chang et al., 2012; Tsai et al., 2012) to detect the specific nucleic acid sequences of *Tritrichomonas foetus*. This reagent set is specially designed to be used on an iiPCR-compatible instrument, **POCKIT™** Nucleic Acid Analyzer. The intended users of the reagent set are veterinarians or technicians who have basic laboratory skills.

This reagent set is intended for research purpose and *in vitro* use only.

SUMMARY AND EXPLANATION

Tritrichomonas foetus (*T. foetus*) is a single celled flagellated protozoan parasite that is known to be a causative agent of bovine trichomoniasis. Bovine trichomoniasis is transmitted through breeding, and it causes the chronic genital tract infection of inflammation and serious problem of reproductive disorder. Generally, bulls infected by *T. foetus* do not show any clinical symptoms (asymptomatic), but can infect cows at mating through directly contact with bull penis or secretions. The parasites initially can infect the vagina, causing vaginitis, and then move to the uterus and oviduct. The serious symptom of cows includes infertility, abortion and delayed or prolonged calving seasons.

PRINCIPLES OF THE PROCEDURE

The assay is based on iiPCR for qualitative detection of *T. foetus*. Fluorogenic probe hydrolysis chemistry is used to generate fluorescent signal when a specific DNA sequence of *T. foetus* is amplified. The primers and probe target 5.8S ribosomal RNA gene of *T. foetus*, and do not cross-react with host genomic DNA and nucleic acids of other pathogens.

PRODUCT DESCRIPTION

A. Materials Provided

- 1) **POCKIT™ *T. foetus* Reagent Set**, 48 tests/package

Component	Contents or Purpose	Amount
Premix Pack	<ul style="list-style-type: none"> ■ <i>T. foetus</i> Premix (lyophilized pellet) containing dNTPs, primers, probe, and enzyme for amplification. ■ Desiccating agent pack. 	6 bags (8 <i>T. foetus</i> Premix vials and 1 desiccating agent/bag)
Premix Buffer B	<ul style="list-style-type: none"> ■ Reaction buffer to re-dissolve the lyophilized pellet. 	2 vials (1.3 ml/vial)
P(+) Control	<ul style="list-style-type: none"> ■ Dried plasmid containing <i>T. foetus</i> partial sequence as positive control. 	1 vial
P(+) Control Buffer	<ul style="list-style-type: none"> ■ Reaction buffer to re-dissolve P(+) Control. 	1 vial (110 µl/vial)
User Manual		1 copy

- 2) **R-tube**, 48 tubes/pack

B. Materials and Equipment Required, but Not Provided

- 1) **PetNAD™** Nucleic Acid Co-prep Kit or **taco™ mini** Automatic Nucleic Acid Extraction System (optional)
- 2) **POCKIT™** Nucleic Acid Analyzer: the iiPCR-compatible instrument.
- 3) **cube™** Mini-Centrifuge
- 4) Micropipette and tips

C. Storage and Stability

- 1) The reagent set should be stored at 4°C and is stable until the expiration date stated on the label.
- 2) Store Premix vials in sealed Premix Pack to avoid hydration of lyophilized components.
- 3) Reconstituted P(+) Control is stable for 6 months at 4°C. Aliquot reconstituted P(+) Control to avoid degradation of nucleic acid.

PRECAUTIONS


- A. Do not open R-tube(s) after reaction to prevent any carryover contamination.
- B. Perform extraction and amplification in two independent spaces to minimize contamination.
- C. Do not reuse R-tube and Premix.

D. Include the P(+) Control to:

- 1) Ensure **POCKIT™** Nucleic Acid Analyzer is working.
- 2) Ensure detection reagent performance after storage.

E. To get optimal fluorescence detection.

- 1) Wear powder-free gloves to handle R-tubes.
- 2) Do not label in the detection area of R-tube.



LIMITATIONS

- A.** The test should only be used for testing nucleic acid extracts. Do not add specimens directly into the Premix.
- B.** **PetNAD™** Nucleic Acid Co-prep Kit or **taco™ mini** Automatic Nucleic Acid Extraction System is recommended for nucleic acid extraction.
- C.** Any deviation from the recommended procedures may lead to sub-optimal results. Performance of the modified protocol should be validated by the users.
- D.** It is strongly recommended to use freshly prepared nucleic acids (within 1 hour after extraction) to achieve optimal results with **POCKIT™** *T. foetus* Reagent Set.

SAMPLE TYPE

This reagent set is intended for analyzing nucleic acids extracted from preputial smegma from bulls and/or vaginal mucus from cows (McMillen and Lew, 2006).

OPERATION PROCEDURE


A. Using POCKIT™ *T. foetus* Reagent Set

- **Note: Before preparing the reactions for iiPCR testing, turn on POCKIT™ Nucleic Acid Analyzer to initiate the calibration for the instrument. The device will complete self-test within 5 minutes. Please refer to the user manual of POCKIT™ Nucleic Acid Analyzer for further details.**
 - **Note: Before using for the first time, add 100 µl P(+) Control Buffer to P(+) Control. Store reconstituted P(+) Control at 4°C.**
- 1) Label R-tube(s) in the labeling area.
 - 2) Prepare one Premix for each sample. (Premix tubes are in Premix Pack. Each Premix Pack contains eight Premix tubes.)
 - **Note: When the pellet is not found at the bottom of the tube, spin tube briefly to bring it down.**
 - 3) Add 50 µl Premix Buffer B to each Premix tube.
 - 4) Add 5 µl nucleic acid extract or dissolved P(+) Control to each

Premix tube. Spin Premix tubes briefly in a mini centrifuge (such as **cube™** Mini-Centrifuge).


- 5) Transfer 50 µl Premix/sample mixture into R-tube.
- 6) Seal top of each R-tube with a cap. Make sure R-tube is capped tightly.
- 7) Place R-tube into the holder of **POCKIT™** Nucleic Acid Analyzer.
- 8) Spin tube/holder set briefly in **cube™** Mini-Centrifuge to make sure all solution is collected at the bottom of R-tube.
 - **Note: Make sure there are no bubbles in the solution.**
 - **Note: Start reaction within 1 hour to reduce the risks of nucleic acid degradation and non-specific reaction.**
- 9) **POCKIT™** Nucleic Acid Analyzer reaction:
 - a) Select "520 nm + 550 nm".
 - b) When "System READY" is displayed, place the holder with R-tube(s) into the reaction chamber.
 - c) Tap cap of each R-tube to make sure the tube is positioned properly.
- 10) Close lid and press "Run" to start reaction program.
- 11) Test results are shown on the monitor after the reaction is complete.

B. Quick Guide



DATA INTERPRETATION

- 520-nm fluorescent signal is used to detect 5.8S ribosomal RNA gene of *T. foetus*; 550-nm fluorescent signal is used for internal control.
- The following example is iiPCR reaction results shown on the monitor.



520 nm	550 nm	Interpretation
+	+	Positive
+	-	Positive
-	+	Negative
-	-	Repeat reaction with freshly prepared nucleic acid.
+	?	Positive
?	+	Repeat reaction with freshly prepared nucleic acid.
-	?	Repeat reaction with freshly prepared nucleic acid.
?	-	Repeat reaction with freshly prepared nucleic acid.
?	?	Repeat reaction with freshly prepared nucleic acid.

ANALYTICAL SENSITIVITY

The detection limit of **POCKIT™** *T. foetus* Reagent Set is about 10 copies/reaction.

TROUBLESHOOTING

Problems	Possible causes	Solutions
False Positive	1) Reuse of micro-centrifuge tubes, tips, R-tubes and Premix.	<ul style="list-style-type: none"> ■ Micro-centrifuge tubes, tips, R-tubes and Premix are for single-use only. Reusing these accessories would cause cross-contamination. ■ Used micro-centrifuge tubes, tips, R-tubes and Premix should be collected and discarded according to local regulation. Do not place the waste close to the working area to prevent cross-contamination.
	2) Contaminated micropipette.	<ul style="list-style-type: none"> ■ Use aerosol-free tips.
	3) Contaminated reagent.	<ul style="list-style-type: none"> ■ Consult with a GeneReach technical support representative or local distributor.
	4) Leakage or spill of reaction from R-tube into reaction chamber of POCKIT™ Nucleic Acid Analyzer.	<ul style="list-style-type: none"> ■ Consult with a GeneReach technical support representative or local distributor.
	5) Contaminated working area.	<ul style="list-style-type: none"> ■ Consult with a GeneReach technical support representative on how to clean up working area.

Problems	Possible causes	Solutions
False Negative	1) Nucleic acid extraction failed.	<ul style="list-style-type: none"> ■ Consult manual of Nucleic Acid Co-prep Set.
	2) Poor nucleic acid quality.	<ul style="list-style-type: none"> ■ Check sample storage condition. ■ Please refer to Troubleshooting section of Nucleic Acid Co-prep Set
	3) PCR inhibition.	<ul style="list-style-type: none"> ■ Do not overload nucleic acid. ■ Spike 5 µl nucleic acid sample into a positive control reaction for a parallel PCR reaction. Negative results indicate the presence of inhibitors in the nucleic acid. In that case, prepare another nucleic acid extract.

REFERENCE

1. Anonymous. Trichomoniasis. In: Merck Veterinary Manual, 9th ed. Whitehouse Station, N.J., Merck and Co. 2005, pp. 1142-3.
2. BonDurant R.H., Anderson M.L. & Blanchard P., (1990). "Prevalence of bovine trichomoniasis among California beef herds". J Am Vet Med Assoc 1990; 196: 1590–1593.
3. BonDurant RH. "Diagnosis, treatment and control of bovine trichomoniasis". Comp Cont Educ Pract Vet 1985; 7: S179–188.
4. Campbell RJ, Schumann F, Janzen E, Wagner B, Polley L. "*Tritrichomonas foetus* in bulls in Saskatchewan". Can Vet J 1993; 34: 571–572.
5. Chang, H. G., Tsai, Y., Tsai, C., Lin, C., Lee, P., Teng, P. A thermally baffled device for highly stabilized convective PCR. *Biotechnology Journal* 2012, 7(5), 662-666. doi: 10.1002/ biot.201100453
6. Felleisen R, Lambelet N, Bachmann P, Nicolet J, Muller N, Gottstein B.. "Detection of *Tritrichomonas foetus* by PCR and DNA enzyme immunoassay based on rRNA gene subunit sequences". J Clin Microbiol 1998; 36: 513–519.
7. Gookin, J.L.; Breitschwerdt, E.B.; Levy, M.G.; Gager, R.B.; Benrud, J.G. "Diarrhea associated with trichomonosis in cats". J Am Vet Med Assoc 1999; 215: 1450-4.
8. Grahn, R.A.; BonDurant, R.H.; van Hoosear, K.A.; Walker, R.L.; Lyons, L.A. "An improved molecular assay for *Tritrichomonas foetus*". Vet Parasit 2005; 127: 39-47.
9. Ho MY, Conrad PA, Conrad PJ, LeFebvre RB, Perez E, BonDurant RH. "Detection of bovine trichomoniasis with a specific DNA probe and PCR amplification system". J Clin Microbiol 1994; 32: 98–104.
10. Kimsey, P.B. Bovine trichomoniasis. In: Morrow, D.A. (Ed.). "Current Therapy in Theriogenology, 2nd ed". Philadelphia, W.B. Saunders 1986, pp. 275-9.
11. McMillen. L., Lew, A.E. "Improved detection of *Tritrichomonas foetus* in bovine diagnostic specimens using a novel probe-based real time PCR assay". Vet Parasit 2006; 141: 204-215.
12. Rae DO, Chenoweth PJ, Genho PC, McIntosh AD, Crosby CE, Moore SA , et al.

- "Prevalence of *Tritrichomonas foetus* in a bull population and effect on production in a large cow-calf enterprise". *J Am Vet Med Assoc* 1999; 214: 1051–1055.
13. Riley DE, Wagner B, Polley L, Krieger JN. "PCR-based study of conserved and variable DNA sequences of *Tritrichomonas foetus* isolates from Saskatchewan, Canada". *J Clin Microbiol* 1995; 33: 1308–1313.
 14. Tsai, Y., Wang, H. T., Chang, H. G., Tsai, C., Lin, C., Teng, P. Development of TaqMan Probe-Based Insulated Isothermal PCR (iiPCR) for Sensitive and Specific On-Site Pathogen Detection. *PLoS ONE* 2012, 7(9), e45278. doi:10.1371/journal.pone.0045278