MegaFLUO® PBS

ad us. vet.
In vitro Diagnostikum

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MegaFLUO® PBS

ad us. vet.
In vitro diagnosticum

PBS (Phosphate buffered saline) tablets for the indirect semiquantitative immunofluorescence detection of specific IgG/IgM antibodies

INSTRUCTIONS FOR USE

6912 Hörbranz – AUSTRIA

1. INFORMATION ON THE TEST-KIT

TEST-KIT COMPONENTS

1 package MegaFLUO® PBS contains:

- 1, 10 or 20 PBS (phosphate buffered saline) tablets, pH 7.2–7.4
- 1 instructions for use

SAMPLE MATERIAL

Serum or plasma

STORAGE AND STABILITY

- The storage temperature for the PBS tablets is +15–25°C
- Shelf life of tablets: 36 months after manufacturing.
- Shelf life of reconstituted solution: 12 months at +4°C

MATERIAL REQUIRED, BUT NOT PROVIDED

MegaFLUO® Slides, Positive Control, Negative Control, FITC anti-species IgG/IgM conjugate, Mounting Medium, washbasin for PBS, test tubes for serum dilutions, microlitre pipettes, 24x50 mm cover slips, fluorescence microscope with filter system for FITC (fluorescein isothiocyanate, excitation wavelength 465–495, barrier filter 515–555) and 400× magnification, 37°C incubator, humid chamber.

INDICATED USE

This product is suitable as reagent for indirect immunofluorescence detection of specific antibodies (IFAT).

LIABILITY

The entire risk due to the performance of this product is assumed by the purchaser. The manufacturer shall not be liable for indirect, special or consequential damages of any kind resulting from the use of this product.

2. TEST PRINCIPLE

Sera are diluted in PBS (pH 7.2–7.4) and dropped onto the slide wells to allow an antigen-antibody reaction at 37°C in the case of a positive sample. By subsequent washing with PBS, non-bound unspecific serum proteins are washed off. In the next step, fluorescein-marked FLUO FITC anti-species IgG/IgM conjugate is added which binds to the antigen-antibody complexes. After an incubation time of 30 minutes, non-bound conjugate is washed off with PBS. Finally, the wells are covered with Mounting Medium and cover slip. Evaluation is done with a fluorescence microscope (filter system for FITC) with 400× magnification.

3. PRECAUTIONS FOR USERS

- Make sure that the test-kit components can be correlated exactly to the particular patient.
- Use a new pipette tip for each sample or dilution.
- The conjugate is photosensitive and sensitive to heat, therefore it should be stored in the dark at 2–8°C.
- The conjugate contains Evans-blue dye, which is potentially carcinogenic. Avoid ingestion and skin contact.
- The sample material, the controls and the slides must be seen as potentially infectious and disposed of accordingly after test procedure, together with the used test-kit components.

4. IMPORTANT ADDITIONAL INFORMATION FOR TEST EVALUATION

Seroprevalence

The cut-off can vary depending on the region and the origin of the sample (dependent on the prevalence and on the state of endemicity). Therefore, it is recommended for every laboratory to determine an individual cut-off.

An acute infection (2- to 4-fold titre increase) can only be determined by the titre determination of a coupled serum test (2 samples in an interval of 2–3 weeks).

The interpretation of test results should always be based on anamnestic and especially clinical data and additional laboratory parameters.

5a. PREPARATION OF THE PBS SOLUTION

1. Dissolve the PBS tablet in ca. 50–100 ml aqua dest by stirring (magnetic stirrer). Adjust the volume to 1000 ml and stir until full dissolution.
2. Sterilise by filtration or autoclaving. Filter the buffer into a sterile flask through a 0.22 µm filter or autoclave for 15–20 minutes and keep at +4°C.
3. The buffer is now ready to use.

5b. TEST PROCEDURE

1. The test-kit components (apart from conjugate!) and the sera to be tested should have room temperature at the time of application.
2. Prepare appropriate dilutions (see instructions for use for the specific test) with PBS for all sera to be tested. For sera found positive on a previous test it is recommended to prepare serial two-fold dilutions in PBS to determine the endpoint titre (= highest dilution that is still positive).
3. Carefully remove the slides from their foil pouch shortly before use and place them into the humid chamber. Apply 1 drop (20 µl) of the Positive and Negative Control on each slide on separate antigen wells. Pipette 20 µl of every serum dilution on separate antigen wells (fig.1). Take care that the antigen wells are completely covered.
4. Incubate for 30 minutes at 37°C.
5. Washing step: Tap remaining serum dilutions gently from the slides and shake the slides gently for 5 minutes in PBS. Repeat this step for another 5 minutes with fresh PBS. Briefly rinse the slides with distilled water. Tap remaining water gently from the slides and, if necessary, dry the teflon mask between the wells with absorbent paper or cotton bud sticks. However, do not allow the antigen wells to dry out!
If using a washing bottle, do not focus the stream directly onto the antigen wells!
6. Place the slides back into the humid chamber and immediately drop 1 drop of FLUO FITC anti-species IgG/IgM conjugate onto each used well. (fig.2). Take care that the antigen wells are completely covered.
7. Incubate for 30 minutes at 37°C and in the dark to protect the photosensitive conjugate.
8. Repeat washing step as described in step 5.
9. Add some drops of Mounting Medium on the cover slips and place them carefully on the slides. Try to remove any possible bubbles carefully.
10. Evaluate the slides using a fluorescence microscope at 400× magnification (fig.3), comparing each well to the fluorescence pattern seen in the Positive and Negative Controls.
11. Sealed slides can be stored at 2–8°C in the dark for up to 7 days.

6. TEST EVALUATION

For the evaluation, a fluorescence microscope with a filter system for FITC (excitation wavelength 465–495, barrier filter 515–555) and 400× magnification is required.

The fluorescence pattern (form, density etc.) of the Negative and Positive Control is considered as reference pattern. Patterns of reactivity different than that seen in the Controls must be considered non-specific, which means negative!

Positive test result

A positive reaction is shown as bright, sharp and clear, yellow-green fluorescence of the pathogen to be detected.

A further dilution of the positive samples is recommended to determine the endpoint titre (= highest dilution that is still positive).

Cut-off fluorescence pattern / recommended cut-off

The pathogen shows a weak (1+) yellow-green fluorescence.

Negative test result

No yellow-green fluorescence visible.

Divergent fluorescent reaction

Reactivity patterns different than the described ones must be considered as non-specific reactions and therefore as negative.

Images are shown at <http://www.megacor.at>

EXPLANATION OF SYMBOLS

- Store at 15–25°C – see label
- For veterinary use only
- In vitro diagnosticum
- Follow instructions for use precisely
- Expiry date – see label
- Lot number
- Do not use test-kit components from different kits, lot numbers or beyond stated expiry date.